

Elektravägen 53 SE-126 30 Hägersten, Sweden +46 10 252 30 00 www.ozonetech.com

# Reliant

The Reliant is an oxygen concentrator ideal for use in small clinics. The Reliant's ability to deliver up to 95.5% oxygen concentration at flows of up to 8 LPM (17 SCFH) make it a great fit for many applications requiring up to 50 psig (345 kPa). An optional module enables the Reliant to supply 8 LPM of compressed, dried, and filtered medical air. Under normal operation, the Reliant is toggled between supplying either oxygen or compressed medical air. The presence of an internal oxygen storage tank and an optional auxiliary oxygen receiver, allows the Reliant to temporarily supply oxygen while simultaneously producing medical air. The duration for such dual mode operation is limited by the size of oxygen storage tanks and the oxygen use rate. This simultaneous dual-supply mode enables the Reliant to briefly provide medical air without an interruption to the flow of oxygen.



#### **Features**

- Self-contained (internal air compressor)
- Microprocessor controlled
- Low operating cost
- Oxygen concentration monitor
- Automatic and unattended operation
- Transportable design
- Easy to install and maintain

### Typical applications

#### Medical

- Anesthesia
- · Oxygen therapy
- Veterinary medicine

#### Available options

- Medical air outlet (factory-installed)
- Auxiliary oxygen receiver (free standing)









## **Specification**

Product characteristics		
Product flow <sup>1</sup>	8 LPM (17 SCFH)	
Product pressure	345 kPa (50 psig)	
Product concentration (nominal)	93%	
Product dew point	-73 °C (-100 °F)	
Dimensions (W x D x H) (nominal)	44x 67 x 91 cm (17 x 27 x 36 in) with flowmeter mounted	
Weight	80 kg (176 lb); 82 kg (181 lb) with optional medical air outlet	
Physical connections		
Product gas outlet	1/4" FNPT x 1/4 NPSM	
Physical connection optional medical air outlet	1/4" FNPT JIC-8	
Ambient operating conditions	Locate the oxygen concentrator system in a well-ventilated area that is protected from weather elements and remains between 4 °C (40 °F) and 44 °C (112 °F)	
Control power requirements (single phase)	120 V $\pm$ 10%, 60 Hz, 6.0 A or 220 V $\pm$ 10%, 50 Hz, 3.0 A	

 $<sup>^1</sup>$  LPM (Liters per minute) gas measured at 1 atmosphere and 21 °C / SCF (Standard cubic foot) gas measured at 1 atmosphere and 70 °F

Elektravägen 53 SE-126 30 Hägersten, Sweden +46 10 252 30 00 www.ozonetech.com







### **Ordering information**

Model	Part number	Description
Reliant	AS132-1	120 V ~ ±10%, 60 Hz
	AS132-2	220 V ~ ±10%, 50 Hz
	AS133-1	120 V $\sim$ ±10%, 60 Hz with medical air outlet [factory-installed]
	AS133-2	220 V $\sim$ ±10%, 50 Hz with medical air outlet (factory-installed)
Accessories		Reliant Standard accessory kit —includes the following:
	KI046-1	• 10 ft stainless braided oxygen hose
		<ul> <li>primary/secondary oxygen ball valve assembly</li> </ul>
	TA150-1	227 Liters/60 gallon oxygen receiver
		Oxygen Reliant Standard accessory kit —includes the following:
	KI048-1	• 7 ft stainless braided oxygen hose
		• oxygen isolation ball valve assembly
		oxygen purge ball valve
Shipping information	Reliant	227 Liter/60 gallon oxygen receiver
Class	92.5	92.5
Commodity classification number	8421.39.8040	7311.00.0000
Dimensions (W x D x H)	71 x 76 x 114 cm (28 x 30 x 45 in) with accessory kit	71 x 76 x 142 cm (28 x 30 x 56 in)
Gross weight	106 kg (233 lb) with accessory kit	80 kg (177 lb)


An unprotected or inadequately ventilated environment, or improper control power may cause damage to the oxygen generator. All performance ratings based on an ambient temperature up to 38 °C (100 °F), up to 304.8 meters (1,000 feet) elevation, and 80% relative humidity.


Elektravägen 53 SE-126 30 Hägersten, Sweden +46 10 252 30 00 www.ozonetech.com











Note: All dimensions are nominal.

Elektravägen 53 SE-126 30 Hägersten, Sweden +46 10 252 30 00 www.ozonetech.com



